

Prevent Falls & Wandering

Quieter CordLess® Fall Monitor without Cords



TL-2100G

Eliminates In-Room Noise

The CordLess® Fall Monitor is cord free and can be mounted where you need it - away from residents - eliminating in-room alarm noise. Monitor can be easily wall-mounted at the bedside or in a hallway up to 30 ft. away using the provided wall-mounting hardware.

Patented Safety Features

This monitor features Patented, cutting-edge technology that automatically lets a caregiver know if the signal between pad and monitor is lost or if the monitor is out of range. In addition, this monitor features an optional Patented TamperProof™ setting, which when enabled requires caregivers to use the provided Caregiver Key™ to reset the alarm. This ensures that residents cannot accidentally reset the monitor.

No Cords!

Cordless alarm system has NO cords between pads and monitors, which reduces tripping hazards and eliminates broken or tangled pad cords.

Key Features:

- **QUIETER** – Fully adjustable volume. Place monitor away from resident (i.e. above bed or in hallway)
- No cords between pads and monitors - reduce tripping hazards and eliminates broken or tangled pad cords
- Patented Pad Signal Lost notification - alerts caregiver if sensor pad and monitor are no longer communicating
- Optional Setting: TamperProof™ Caregiver Key allows control of who is able to reset the monitor
- Raised Visual Status light enabling caregiver to see alarm from a distance
- Monitor up to two components at one time
- Batteries not included (3-AA)
- Optional AC adapter available (AC-05)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com