

# Prevent Falls & Wandering

## Low Cost, Quiet and CordLess Fall Prevention Alarm



**No Cords! Place monitor where you need it!**

### Low Cost Solution

The Economy CordLess® Exit Alarm provides affordable and effective cord-free monitoring with a range up to 300 ft. Ideal for home or small facility use.

### Easy To Use

The Economy CordLess® Exit Alarm is easy-to-use and works with all CordLess® Sensor Pads & Floor Mats, Nurse Call Buttons, and Motion Sensors. Monitor is light-weight and can be paired with up to six (6) wireless components to serve as a portable caregiver alert.

### Remove Noise From the Room

Monitor can be mounted away from the bed, removing in-room alarm noise and reducing entanglements and tripping hazards. Adjustable volume and gentle chime alert won't startle residents and offers QUIETER fall prevention.



433-EC

### Key Features:

- Affordable & effective wireless monitoring
- Easy set up and operation
- Monitor up to six (6) components with one alarm
- Full Range Volume Control
- Batteries not included (3-C)
- Optional AC adapter available (AC-04)



**Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)**