

# Prevent Falls & Wandering

## Reduce Risk of Falls, While Increasing Caregiver Freedom

Use quiet monitor with wireless signal to pager



### Reduce Alarm Noise in Room!

The monitor can be set to silent to reduce in-room alarm noise. When a resident attempts to get up, a silent wireless signal is sent to a caregiver pager or central nurse's station.

### Flexible & Versatile

The wireless fall monitor is designed to work with Corded Bed & Chair Sensor Pads, Floor Mats, Seat Belts, Pagers, Emergency Call Lights, Nurse Call Cords and Central Monitoring Units. This monitor allows caregivers the flexibility to monitor residents based on their unique needs.

### Increase Caregiver Freedom!

Caregiver can attend to other duties and receive an alert with use of a pager when someone is getting up from a bed or chair. Can also be used with a floor mat to notify if they are approaching a door.

TL-2016  
Quiet Wireless  
Fall Monitor

TL-2016P  
Caregiver Pager

### Key Features:

- Monitor sends wireless signal to pager up to 300 feet away
- Call Button on front of monitor
- Nurse Call Capability - enables you to use with existing hard wired nurse call system
- Batteries not included (3-AA)
- Optional AC adapter available (AC-05)



Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)