

Prevent Falls & Wandering

Eliminate The Need to Check Dates on Sensor Pads



TL-2100CP

Keep Track Of Dates On Pads

The TL-2100CP automatically tracks pad life so caregivers no longer have to write dates on pads. The Monitor counts down the days the pad has been in use and the Check Pad Light will illuminate 10 days prior to pad expiration. A short beep will also alert caregivers that it's time to change pad. When used with Timed Bed & Chair pressure pads, this monitor will not only save time and staff resources, but also **Eliminate F-Tags!**

Easy-To-Use

Monitor can be placed at the bedside or on a wheelchair and connected to a corded sensor pad. The sensor pad is placed under the resident and when the resident gets up and pressure is removed from the pad, the monitor will sound alerting the caregiver. Monitor will automatically reset when pressure is re-applied to the connected sensor pad.

Patented TamperProof™ Technology

This monitor has been designed with an optional TamperProof Setting which, when enabled, requires caregivers to use the provided Caregiver Key to reset the alarm. This prevents resident tampering and ensures that only the designated caregiver can silence alerts.

***Must be used with Timed Sensor Pads**

Key Features:

- Full-Range Volume Control
- Nurse Call Capability - enables you to use with existing hard wired nurse call system
- Optional Setting: TamperProof™ Caregiver Key allows control of who is able to reset the monitor hard wired nurse call system
- Batteries not included (3-AA)
- Optional AC adapter available (AC-05)



Call Us Today! (800) 650-3637 | www.smartcaregiver.com