

# Prevent Falls & Wandering

## Fall Alarm That Works With Bed Pad & Pull-String



TL-3100V

### Unique Dual Purpose Fall Monitor

The TL-3100V DUAL Recordable Voice Alarm works as a pull-string monitor OR sensor pad monitor and can easily be transferred from bed to wheelchair and back. This eliminates the need for multiple monitoring devices. The alarm alerts when the magnetic pull-string is pulled free of the alarm or when pressure is removed from weight-sensing pad. Alarm resets when the magnetic pull cord is replaced, when pressure is re-applied to the connected sensor pad or when reset button is used.

### Personalized Recordings

Eliminate startling alarm noise by recording a personalized message to be played back when the alarm is activated. With up to 15 seconds of playback time available, you can record a melody or familiar voice to be heard as the alert sounds.

### Key Features:

- Professional Grade
- Dual Functions
- Recordable Voice Option
- Pad Lost, Low Battery & Alarm Status Indicator Lights
- Nurse Call Capability
- Monitoring Options: Corded Bed Pad, Chair Pad, Floor Mat, Seat Belt, and Fall Mattress
- Optional AC adapter available (AC-05)
- Batteries not included (3-AA)



Call Us Today! (800) 650-3637 | [www.smartcaregiver.com](http://www.smartcaregiver.com)