

Prevent Falls & Wandering

Low Cost Window/Door Exit Alarm

Sends Signal To Central Monitor or CordLess® Alarm



433-CMU
Economy Central
Monitoring Unit



433-EC
Economy CordLess®
Alarm

Know when a window or door has been opened

Easy To Use

The 433-EXT Window/Door Exit Alarm can be mounted to a window or door frame using the provided mounting hardware or adhesive tape.

Alerts Caregiver

When the window or door opens the contact magnet is separated from the sensor, a wireless signal is sent to the Economy Central Monitoring Unit or Economy CordLess® Alarm to alert the caregiver.

Program Multiple Exit Alarms

Multiple Door/Window Exit Alarms may be programmed to the Economy Central Monitoring Unit or Economy CordLess® Alarm to monitor more than one area of a facility or home at a time.

Battery included (1-12V)



433-EXT
Low Cost Wireless
Window/Door Exit Alarm



Call Us Today! (800) 650-3637 | www.smartcaregiver.com